



UMS *Sculptor*

Newest Second Generation UMS up to 150HZ

Muscle Building

Lipo-Reduction

Body Shaping

Vaginal Relaxation

Pelvic Floor Muscle Relaxation





Patented handles with two types: flat and curve designed 2 or 4 handles can work at the same time, meeting the simultaneous operation to different clients and different target areas

UMS *Sculptor*

REDEFINING YOUR REAL SELF

Muscle Building

Pulsed Electro Magnetic Wave

High technology exclusively dedicated to a slimmer body



NON-INVASIVE · SAFE · EFFECTIVE



OPERATOR INDEPENDENT



IMORE COSMETICS · ENHANCE RESULTS




UMS Sculptor


Equipped with 4 worldwide unique applicators, 2 different sizes worldwide

Unique small handle curved, ergonomic, perfect for arm, waist, hips, thigh and calf.



UNIQUE SMALL HANDLE CURVED





Muscle anatomy

HOW IT WORKS ?

IT IS NOT ONLY ABOUT FAT



Voluntary contractions may increase the demand for energy from fat cells. Released epinephrine signals the fat cells to initiate lipolysis.



Intensive supramaximal contractions enhance release of epinephrine which triggers a cascade effect leading to supramaximal lipolysis in fat cells.




The free fatty acids (FFAs) overflow causes cell dysfunction and induction of apoptosis - programmed cell death. Dead cells collapse and are naturally flushed away.




Intensive supramaximal contractions stimulate muscles' response with a deep remodeling of its inner structure, results in increased muscle density and volume.






Myofibrillars loosen

As age grows, Myofibrillars are loose and fragile because of pregnancy, fat accumulation, lack of exercises and other factors. People become fatter and flabbier.



Myofibrillars break

The equipment interact with muscle neurons to induce supra maximal muscle contractions. In the process of continuous muscle contractions, part of Myofibrillar are gradually broken to cause muscle self-repair.




Myofibrillars repair

After treatment, Myofibrillars tend to be strong and tight. Muscle becomes firm and elastic. It responds with a deep remodeling of inner structure that result in muscle building and growth and body shaping.

TO GET AN ATTRACTIVE BODY YOU NEED TO BUILD MUSCLE

UMS Sculptor



treatment, 30 minutes,


20000 crunches or squats

100% supramaximal muscle contractions,

other exercises only reach to 20-30%


Abdominal muscle mass increased by 16%

Abdominal muscle firmness increased by 11%




Before treatment:

Muscle is thin and the circumference of muscle is short. Myofibrillar is loose with small density.




After treatment:




Muscle is strong and the circumference of muscle is long. Myofibrillar is strengthened and increased with large density.



Myofibrillar



Myofibrillar increased firmly



FAT ACCUMULATION

When gain excess weight, more fat accumulates in the existing fat cells, causing them to swell and form cellulite. Fat accumulation can occur in the abdomen, buttocks and thighs. become fatter and flabbier.

FAT DISRUPTION


Electro Magnetic Wave, intensive supramaximal contractions enhance release of epinephrine which triggers a cascade effect leading to supramaximal lipolysis in fat cells.

FAT LOSS

After UMS Sculptor treatment, targeted fat cells shrink in volume and are eliminated.


Before treatment:

Subcutaneous fat and visceral fat are thickness. Muscle is thin with fatty and loosen body.




After treatment:

Subcutaneous fat and visceral fat are reduced. Muscle is strengthened and increased with firm and strong body.



TO GET AN ATTRACTIVE BODY YOU NEED TO BURN FAT

UMS Sculptor



Beautiful+ Health

Treatment/ Week

One course of treatment, Reduction of fat is 19%


Waist circumference loss 4cm

Fat metabolism increased to 5times



UMS Sculptor


NON-INVASIVE BUTTLIFT




MUSCLE TONING
NON-INVASIVE BUTTLIFT




BUTTOCKS PROBLEMS




M Type




After
Childbirth



Wide Hips



V Type



A Type


Ideal Type
Peach Buttocks

Effect



BUTT, OF COURSE....

Introducing **UMS Sculptor**
a non-invasive treatment, is non-ionizing,
non-radiating and non-thermal.
Best of all,
UMS Sculptor is anesthesia
and needle free.





UMS Sculptor


ABDOMEN SCULPTING




FAT BURNING
ABDOMEN SCULPTING



ABDOMEN PROBLEMS



Effect



BEFORE


Skin, fat and muscle form your overall body appearance.

DURING

The skin remains unaffected while the energy penetrates to fat and muscle layers.

AFTER

The skin remains unaffected while the energy penetrates to fat and muscle layers.



UMS *Sculptor*

ARMS AND THIGHS FIRMING



Arms Firming Solution




Thighs Firming Solution

Firm the skin:

Due to weight fluctuations, pregnancies, passing time, and other factors, the skin gradually loses its tonicity and suppleness. Even though this skin sagging affects the entire body, certain areas are more susceptible to it: arms and inner thighs.

UMS Sculptor

PRIVATE HEALTH OF WOMEN




Rectus Abdominis Firming Solution

Firm the skin:

UMS Sculptor uses high-energy focused electromagnetic wave technology and high-magnetic pulses to generate directional motion force. One operation is equivalent to 10,000 Kegel exercises. It stimulates motor neurons to control pelvic floor muscle contraction, activates pelvic floor tissue, and promotes collagen regeneration.


Indications: vaginal relaxation, pelvic floor muscle relaxation, urinary incontinence, uterine prolapse, etc.




Vagina Firming Solution

MUSCLE BUILDING
AND FAT BURNING
NO SWEAT REQUIRED!

UMS Sculptor



Traditional
Muscle
Building



Other Body
Shaping
Devices



Surgical
Liposuction

Specification UMS Sculptor TECHNICAL DATA

Technology: High Energy Pulsed Electromagnetic
Applicators: 4 Applicators
Magnetic Power: 12.0 Tesla \pm 20%
Magnetic Frequency: Symmetric Biphasic Pulse, Max
150Hz

Treatment Time: 30 Minutes
Pulse Width: 260 μ s
Interface: 10" Touch Screen
Power: 3500W
Dimension: 460 \times 380 \times 1130 (mm)

Accessory

